

# FlowerGraft™ Preshaped Bone Wedge System

FlowerGraft™ Preshaped Bone Wedge System is a complete line of preshaped structural allografts for Evans and Cotton procedures. FlowerGraft Preshaped Bone Wedges offer surgeons convenient, safe, reliable, Ready-for-Surgery™ options for performing reconstruction procedures in the foot and ankle.

## Surgically Efficient

Allograft bone wedges are provided sterile, stored at ambient temperature, packaged in saline and do not require hydration prior to use.

The Ready-for-Surgery™ procedure kits contain all the instruments necessary to select & insert the appropriate graft

Anatomic preshaped wedges do not require intra-operative shaping.

## Smart Design

Bicortical faces provide structural support and cancellous core provides an excellent scaffold for bone growth.


## Plating Options

Low profile, procedure specific plates provide unique anatomic fit and accommodate the range of available graft options.


# FlowerGraft™ Preshaped Evans & Cotton Wedges

The Evans and Cotton Wedges are available in the following sizes to accommodate anatomic variability.

| Part #  | Description                         | H    | H1  | L    | W    |
|---------|-------------------------------------|------|-----|------|------|
| AGB 106 | FlowerGraft™ Pre-Shaped Evans Wedge | 6mm  | 4mm | 22mm | 22mm |
| AGB 108 | FlowerGraft™ Pre-Shaped Evans Wedge | 8mm  | 5mm | 22mm | 22mm |
| AGB 110 | FlowerGraft™ Pre-Shaped Evans Wedge | 10mm | 5mm | 22mm | 22mm |
| AGB 112 | FlowerGraft™ Pre-Shaped Evans Wedge | 12mm | 5mm | 22mm | 22mm |



| Part #  | Description                          | H   | L    | W    |
|---------|--------------------------------------|-----|------|------|
| AGB 305 | FlowerGraft™ Pre-Shaped Cotton Wedge | 5mm | 22mm | 16mm |
| AGB 306 | FlowerGraft™ Pre-Shaped Cotton Wedge | 6mm | 22mm | 16mm |
| AGB 307 | FlowerGraft™ Pre-Shaped Cotton Wedge | 7mm | 22mm | 16mm |
| AGB 308 | FlowerGraft™ Pre-Shaped Cotton Wedge | 8mm | 22mm | 16mm |



# Evans & Cotton Osteotomy Procedure Kits

The Evans Osteotomy Procedure Kit (AGB 901) and Cotton Osteotomy Procedure Kit (AGB 902) are provided sterile and consist of

- ▶ Trial Sizer – For selecting appropriate sized graft
- ▶ Tamp – For inserting graft into osteotomy site
- ▶ K-Wires (1.6mm dia) – For provisional fixation

## Evans Osteotomy Procedure Kit



## Cotton Osteotomy Procedure Kit

◀ TRIAL SIZER ▶




◀ K-WIRES ▶

◀ TAMP ▶

## FlowerADC™ Kit (optionally available instrument)


The FlowerADC™ Kit maybe used to Distract/Compress and maintain stability in osteotomy procedures.



## Evans Plates (optional)

3-Hole plates are provided sterile and are available in multiple sizes to account for anatomic variability.

| Part #  | Description               | W    |
|---------|---------------------------|------|
| AGB 116 | Evans 3-Hole Plate, Right | 16mm |
| AGB 120 | Evans 3-Hole Plate, Right | 20mm |
| AGB 216 | Evans 3-Hole Plate, Left  | 16mm |
| AGB 220 | Evans 3-Hole Plate, Left  | 20mm |



©2017 Flower Orthopedics. All rights reserved. FlowerCube, FlowerCarriage, FlowerDerm, FlowerGraft, IndicationCube, IndicationCarriage, FlowerADC, FlowerCart, Ready-for-Surgery, FlowerPatch, FlowerFlo and FlowerPot are trademarks of Flower Orthopedics. Please refer to package insert for description, indications, contraindications, warnings, precautions, and other important information.



Flower Orthopedics Corporation • 100 Witmer Road, Suite 280 • Horsham, PA 19044  
T: 877.778.8587 • F: 215.394.8904 • [www.flowerortho.com](http://www.flowerortho.com)

FOC-FLWRGRFT 05/17